



Food-Mouth-Disease Virus Full line-up system


LiliF® Diagnostic System provides a “ Total solution from Screening to Confirmation ”

More Simple, More Accurate, More Effective

● LiliF® Full line-up system in FMDV Diagnosis

Rapid Screening System

Easy portable
User-friendly
Decreased turnaround time



Extraction & Automatic System

Stable extraction of DNA/RNA from
pathogenic samples
Spin column & Automatic type



Molecular Diagnostics

High Specificity & Sensitivity
Patent Technologies



● LiliF® FMDV Products Information



❖ LiliF® FMDV Ab rapid test kit

- ✓ LiliF® FMDV Ab Rapid kit is a Lateral Flow Immuno-Chromatographic Assay kit for the rapid and accurate qualitative detection of antibody against FMDV in the blood.



❖ LiliF® FMDV Real-time RT-PCR Kit

- ✓ LiliF® FMDV Real-time RT-PCR Kit is able to detect directly and specifically FMDV polyprotein by Real-time PCR technology on the basis of a genetic database of target nucleic acid fragments

● Related Products

Type	CONTENTS	CAT NO	Size	Specimen / Target Gene
Ab test	LiliF® FMDV Ab test	IR0131176	10 tests	Wh, S, P
Real-time RT-PCR	LiliF® FMDV Real-time RT-PCR Kit	IPC11031	50 tests	VP-1 (Target gene)