


# Food Testing Solutions by Real-time PCR

**3 types** of  
Food Testing Solutions

Washington, D.C., June 7, 2021 - "**Safe food today for a healthy tomorrow**" is the theme of this year's World Food Safety Day, which promotes actions that help prevent, detect and manage foodborne risks.

(Credit : WHO)

## 2021 Recalls by Category




(Credit : US CDC)







**LiliF® Salmonella Real-time PCR Kit** is designed as a probe-based Real-time PCR product targeting the specific sequence in the *invA* gene of Salmonella.



**LiliF® Listeria Real-time PCR Kit** is designed as a probe-based Real-time PCR product targeting the specific sequence in the *hly* gene of *Listeria monocytogenes*.




**LiliF® VTEC Real-time PCR Kit** is able to detect directly and specifically *E. coli* all serotype(H1 to H18) by real-time PCR technology on the basis of a genetic database of target nucleic acid fragments.



## Characteristics

Hot Start Function / High Sensitivity & specificity / Fast & Accurate

## Work Flow



## Order Information

DETECTION	Product	Cat No.	Capacity
	LiliF® Listeria Real-time PCR Kit	IP24386	50 Tests
	LiliF® Salmonella Real-time PCR Kit	IP24632	50 Tests
	LiliF® VTEC Real-time PCR Kit	IP24736	48 Tests
EXTRACTION	Product	Cat No.	Capacity
	Miracle-AutoXT Automated Device	IMC-NC15PLUS	1 Unit
	AutoXT total gDNA kit	17188-96/48	96/48 Tests
	G-spin total DNA Extraction kit	17045	50/200 col

• If you have any questions, please send an email to us. [tech@intronbio.com](mailto:tech@intronbio.com)